









Multi-Unit | Straight & Angulated

The Multi-Unit system provides a solution for screw-retained prostheses even on complicated-to-restore implants (for example, multiple tilted implants). The Multi-Unit system comprises a full range of sizes for both the upper and lower jaws. Straight, 17°, 30° and 45° adaptors, in a variety of heights, connects to a wide range of complementary products.

Material: Titanium (Ti6Al4V ELI)



Straight	L (mm)
 NM-A7101	1
 NM-A7102	2
 NM-A7103	3
 NM-X7100	Included with all multi-unit bases

Angulated 17°	Angulated 30°	Angulated 45°	L (mm)
 NM-A7112	-	-	2
 NM-A7113	 NM-A7133	-	3
-	 NM-A7134	 NM-A7144	4
 NM-S7101	 NM-X7101	Included with all multi-unit bases	

Recommendation: Tighten the base at a torque of 25 Ncm.

Recommendation: Tighten the screw at a torque of 20 Ncm.