









Vari-Connect

The Vari-Connect system presents a complete solution for removable prostheses on tilted implants. It provides all the required equipment for removable prostheses, both on ball attachments and flat attachments, covering a wide range of possible situations. Straight, 17°, 30° and 45° adaptors are available. Complementary products are fixed to the adaptors by the adaptor thread.

Material: Titanium (Ti6Al4V ELI)



Angulated 17°	Angulated 30°	Angulated 45°	L (mm)
 NM-A7212	-	-	2
 NM-A7213	 NM-A7233	-	3
-	 NM-A7234	 NM-A7244	4
 NM-S7101	 NM-X7101	Included with all vari-connect bases	

Recommendation: Tighten the screw at a torque of 30 Ncm.

Vari-Connect Components

Healing Cap



Screw
NM-S7102

Healing Cap
Ø 4.8mm, H 4.3mm
NM-H7101

Transfers & Analog



Screw
L 13mm
NM-S7111

Open Tray
Ø 4.8mm, H 11mm
NM-T7111



Plastic Cap
Ø 4.8mm, H 11mm
NM-T4402


Screw
NM-S7102

Closed Tray
Ø 4.8mm, H 7mm
NM-T7101


Closed Tray
Ø 4.8mm, H 8mm
NM-T7102

Analog
Ø 4.8mm
NM-T7251


Ball-Connect




Ball-Connect
Ø 4.8mm, H 2mm
NM-T7262




Ball-Connect
Ø 4.8mm, H 3mm
NM-T7263




Ball-Connect
Ø 4.8mm, H 4mm
NM-T7264




Metal Cap
Ø 5mm, H 3.2mm
NM-T3001



Silicon Extra Soft Cap
NM-T3002




Silicon Soft Cap
NM-T3003



Silicon Standard Cap
NM-T3004


Star Hex. Drivers



Star Hex. Driver
☆ 1.25mm, L 7mm
NM-X7006



Star Hex. Driver
☆ 1.25mm, L 14mm
NM-X7007



Long Driver
☆ 1.25mm, L 20mm
NM-X7011

Flat-Connect



Flat-Connect
Ø 4.8mm, H 3mm
NM-T7273



Flat-Connect
Ø 4.8mm, H 4mm
NM-T7274



Metal Cap
Ø 5.5mm, H 2.3mm
NM-T3010



Silicon Extra Soft Cap
NM-T3015



Silicon Soft Cap
NM-T3016



Silicon Standard Cap
NM-T3017




Laboratory Cap
NM-T3018




Sealing Disc
NM-T1824


Star Hex. Motor Mounts



Motor Mount
☆ 1.25mm, L 20mm
NM-X7008



Motor Mount
☆ 1.25mm, L 25mm
NM-X7009



Motor Mount
☆ 1.25mm, L 35mm
NM-X7010